JUN-AIR®

Compressor - Vacuum pump

**Model OF301 DC motor / OF301 motor / OF311 motor / OF301 vacuum
OF302 motor / OF312 motor / OF322 motor / OF302 vacuum
OF301-4B / OF302-4B / OF302-15B / OF302-25B / OF302-25BD2
2xOF302-40B / 2xOF302-40BD2**

Operating manual

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Mode d'emploi

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Betjeningsforskrift

JUN-AIR®

GB	Operating manual	5
DE	Betriebsanweisung.....	10
FR	Mode d'emploi.....	15
ES	Modo de empleo	21
NL	Gebruiksaanwijzing.....	26
DK	Betjeningsforskrift	31
	Technical data	37
	Technische Daten	37
	Caracteristiques techniques	37
	Detalles técnicos.....	37
	Technische gegevens	37
	Tekniske data	37
	Diagrams	47
	Zeichnungen	47
	Dessins ..	47
	Diagramas	47
	Tekeningen.....	47
	Diagrammer	47
	Spare parts.....	53
	Ersatzteile	53
	Pièces détachées	53
	Piezas de recambio	53
	Onderdelenlijst.....	53
	Reservedelsliste	53
	Pictures/illustrations.....	65
	Abbildungen/Illustrationen	65
	Photos/illustrations.....	65
	Fotos/ilustraciones.....	65
	Afbeeldingen/illustratien	65
	Billeder/illustrationer	65

Operating manual

Information

Please note that you can find the pictures and figures we are referring to on page 65.

Important - read this first!

Please read the following information and operating instructions included with this product before use. This information is for your safety and it is important that you follow these instructions. It will also help prevent damage to the product. Failure to operate the unit in accordance with the instructions or using JUN-AIR unauthorized spare parts can cause damage to the unit and could cause serious injury.

CAUTION: To reduce risk of electric shock

- Only authorized service agents should carry out service. Removing parts or attempting repairs can create an electric shock. Refer all servicing to qualified service agents.
- If this unit is supplied with a three-pin plug, connect with a properly earthed outlet only.

WARNING: To reduce risk of electrocution

- Do not use this unit with electrical voltages other than stated on the rating plate.
- Never leave this product unattended when plugged in.
- Always unplug this unit immediately after use.
- Store in a dry place.
- Do not use this product in or near liquid or where it can fall or be pulled into water or other liquids.
- Do not reach for this product if it has fallen into liquid. Unplug immediately.
- This unit is not weatherproof. Never operate outdoors in the rain or in a wet area.

DANGER: To reduce risk of explosion or fire

- During spraying with combustible liquids risk of explosion may arise, particularly in closed rooms
- Do not use this product in or near explosive atmospheres or where aerosol products are being used.
- Do not pump any other gases other than atmospheric air.
- Do not pump combustible liquids or vapours with this product; do not use it in or near areas with combustible or explosive liquids or vapours.
- Do not use this unit near naked flames.

CAUTION: To prevent injury

- Compressed air can be dangerous; do not direct airflow at a persons head or body.
- Always keep the compressor out of reach of children.
- Never operate this product if it has a damaged power lead or plug, if it has been dropped or damaged, or if it has fallen into water. Return the product to a service centre for examination and repair.
- Keep the electrical cable away from hot surfaces.
- Ensure all openings are kept free of restriction and never place the motor on a soft surface where the openings may be blocked. Keep all openings free from dust, dirt and other particles.
- Never insert fingers or any other objects into fans.
- This unit is thermally protected and can automatically restart when the overload resets.
- Wear safety glasses, when servicing this product.
- Use only in well ventilated areas.
- This product may only be connected to units or tools with a max.

pressure higher or equal to that of the compressor.


- The surface of the compressor can get hot. Do not touch compressor motor during operation.

Failure to observe the above safety precautions could result in severe bodily injury, including death in extreme cases.

IMPORTANT: General directions for use

- Protect compressor against rain, moisture, frost and dust.
- The compressor is constructed and approved for a max. pressure as stated under Technical Specifications.
- Do not operate compressor at ambient temperatures exceeding 40°C/104°F or falling below 0°C/32°F.
- If the supply lead on the compressor is defective, an authorized JUN-AIR distributor or other qualified personnel must carry out the repair.

Warranty

Provided that the operational instructions have been carried out, your JUN-AIR compressor is guaranteed against faulty material or workmanship for 2 years.

The air receiver is guaranteed for 5 years.

The guarantee does not cover damage caused by violence, misuse, incorrect repairs or use of unoriginal spare parts.

Costs of transportation of parts/equipment are not covered by the guarantee.

JUN-AIR's Conditions for Sale and Delivery will generally apply. JUN-AIR International A/S reserves the right to change technical specifications/constructions.

About this manual

The instructions in this manual are intended for the OF300 compressors and vacuum pumps and most of the text applies to both options. However, when this is not the case the following symbols are used:

= applies to compressor only

= applies to vacuum pump only

Contents of box (motor only)

Your new compressor/vacuum pump should be delivered in a clean and undamaged box. If not, contact your distributor immediately. The box should contain the following:

- 1 OF300 motor
- Kit
- 1 capacitor
- 1 operating manual

Installation (motor only)

Your JUN-AIR motor is very easy to operate. Observe the following simple instructions and you will get many years' service from your motor.

- Visually inspect unit for shipping damage, contact your supplier immediately if you think the unit may have been damaged.

Warning!

To avoid the risk of electric shock, do not install the compressor in areas where it may get in contact with water or other liquids. In which case protection is required.

- Remove the plastic protection from the outlets. Mount the enclosed O-rings on the plugs supplied with the motor in the outlets that are not used (fig. 1a). Mount the nipple in the outlet port, secure with loctite 275 and connect the pressure pipe (fig. 1b + fig. 1c). Use only pressure pipes tolerating a constant temperature of min. 150°C. Ensure that the pressure pipes have a sufficient internal diameter to avoid pressure loss in the system. Check for leaks.
- Mount the supplied feet on the compressor (fig. 2). It is possible to change the distance between the feet depending on the actual mounting of the compressor (see drawings with various feet measurements). It is important that the compressor is mounted vertically to avoid damage to the rubber feet. If the compressor is to be mounted in another position, please contact JUN-AIR International A/S for further information. The mounting of the rubber feet in the 90 mm track will reduce vibrations.
- Sufficient cooling from the surroundings is important. Place the motor in a dustfree, dry and cool, yet frostfree, room. Do not install in a closed cupboard, unless adequate openings for ventilation are available on top and bottom (minimum 500 cm²/77.5 in² each). If the motor is placed under a table, a minimum of 10 cm/3 inch free height must be available above the motor (fig. 3), or an opening of Ø30cm/11.8 inches, corresponding to the top of the motor, may be cut in the table. Ensure that the motor stands firmly on the floor.
- The intake air may be supplied from another place (for instance outside).

Nipple and hose should then be mounted in the thread on the top of the intake filter, and led to the external air source. Ensure that hoses for the intake are sufficient to avoid performance loss.

Electrical installation (motor only)

Warning!

Incorrect electrical connection may result in electric shock. The electrical connection must be carried out in accordance with local electrical regulations and by qualified electrical engineers.

Note!

Earthing of all AC models must be ensured during installation. The capacitor must be earthed, as failure to do so may cause electric shock when touched. Plug the motor into an earthed socket of nominal voltage and ensure that fusing is adequate, see Technical Specifications page 35.

AC models

- For electrical connection, refer to page 44.
- Check motor plate for frequency, voltage and capacitor and ensure that it corresponds to the voltage and frequency used for the motor. The voltage stated on the motor plate: 120/240V (/) means that the motor can operate at 120V or 240V, but this requires recoupling of the internal wiring from the electrical motor (see the electrical diagrams). 220-230V (-) means that the motor may operate within the range of 220V to 230V without recoupling of the internal wiring.
- Mount bracket underneath the motor (fig. 4a).
- Mount capacitor vertically on mounting bracket (fig. 4b).

DC models

- For electrical connection, refer to page 45.
- During installation of the DC motor, check the motor plate for the direction of rotation as it is important that the "+" and "-" poles are connected in a correct way (fig. 5). Otherwise, the fan may rotate in the wrong direction and the motor will not be cooled. Consequently, the fan will not cool the motor efficiently, causing a shorter lifetime of the compressor, and invalidating warranty.
- A DC motor may be coupled directly on to the "+" and "-" poles without thermal protection. In case of high temperatures and leaks, we recommend the use of an external relay. This relay will switch off the motor if the motor is thermally overloaded (see electrical diagram for DC motors).

Operation (motor only)

- If the temperature of the motor is extremely low (for instance after transportation or stocking), allow motor to get to room temperature before switching on the motor.
- Do not use motor for compression of liquids and dangerous gasses, such as petrol vapour and solvents.

Important!

The compressor/vacuum pump is only suitable for atmospheric air.

- Do not remove protection covers during operation as it may cause electric shock or risk of other personal injury (fig. 6).
- Ensure that motor is correct for air supply/vacuum flow required, see Technical Specifications, page 37.
- If supplied with a receiver, open the outlet cock on the receiver and connect equipment.

Preparing the motor (motor only)

Mounting of vibration dampers OF301 - 50Hz

- Mount the unbraco bolts and nuts in holes on dampers, as shown (fig. 7).
- Tighten the bolts when the dampers have been correctly placed, by the protection covers at both ends of the motor as shown (fig. 8).
- Fasten the motor to the correct mounting, using 4 M6*10 bolts and nuts.

Adjustment of pressure switch (comp.)

Warning!

If maximum pressure is exceeded, reduced lifetime may result. Contact JUN-AIR for information on operation at higher pressure.

- All AC and DC motors may run at 100% continuous operation, but 50% operation is recommendable to prolong lifetime.
- Do not lubricate the oil-less motor with oil, as this will destroy important components.
- OF302 may be delivered with one cylinder for vacuum and one for compression, please contact JUN-AIR for further information.
- Start the compressor using the 0/1 switch on the pressure switch (fig. 9). The compressor will automatically switch off at the preset pressure. If the motor does not start it may be due to pressure in the receiver, and the motor will then start automatically when the pressure reduces to approx. 6 bar/87 psi.
- Adjustment of pressure (fig. 10):
A: Max. pressure adjustment (cut-out)
B: Differential adjustment (cut-in)
The cut-in pressure (normally 6 bar) is set by adjustment of differential screw B. Turn clockwise to reduce cut-in pressure. The cut-out pressure is set by even adjustment of the two screws A. (Cut-in pressure + differential = cut-out pressure). Turn clockwise to increase cut-out pressure. The switch is normally factory set for operation at 6-8 bar (approx. 90-120 psi).

Adjustment of vacuum switch (vacuum)

Note!

Vacuum switch is only mounted on units with receiver.

Before adjusting the vacuum switch, adjust the left screw, see (fig. 11) (turn clockwise), so that the indicator on the left scale is all the way to the top and adjust the right screw, so that the indicator on the right scale is all the way at the bottom as shown (fig. 12).

This base adjustment results in continuous operation and a minimum differential pressure of the vacuum switch, e.g. approx. 150 mbar.

Subsequently, the exact adjustment may take place.

Adjustment of maximum vacuum (left screw)

Start vacuum - Differential vacuum = stop vacuum

The procedure below is based on the base adjustment

described above.

1. In practice, the vacuum pump is to operate until the vacuum required is reached (see vacuum gauge).
2. Switch off the vacuum pump by turning the left screw counter clockwise.
3. Open the air inlet on the receiver in order to ensure automatic start of the vacuum pump.
4. Close the air inlet in order to ensure automatic stop of the vacuum pump at vacuum required.
5. Continue the process above (1-4) until the vacuum pump stops at vacuum required.

Example:

If the vacuum required is e.g. 300 mbar, the vacuum pump will not start operating until the vacuum is approx. 450* mbar. (The difference of 150 mbar corresponds to the adjustment on the differential scale).

* 450 mbar vacuum corresponds to -0.45 on the vacuum gauge.

Fault finding and repair (all)

Important!

Switch off and isolate from electrical supply before removing any parts from the pump. Empty air receiver of air before dismantling parts of compressor unit's pressure system.

1. Compressor does not start:

- a) No power from mains. Check fuses and plug.
- b) Breakage or loose joints in electrical connections.
- c) Defective capacitor.
- d) The thermal protection has switched off the pump due to overheating. When cooled the pump will automatically turn on at a suitable operation temperature. Go through the points in section 6.
- e) The compressor has not been unloaded and there is back pressure on the piston. Make sure that the compressor is unloaded each time it stops.
- f) The pump is locked.
- g) Pressure in the air receiver is too high for activation of the pressure switch. The pressure switch makes circuit only when pressure has dropped to preset start pressure. Empty receiver.

⚠ If the vacuum pump is fitted with receiver and vacuum switch:

- h) Vacuum in the receiver. Vacuum pump will switch on when vacuum reaches the minimum vacuum preset on the vacuum switch. Open the inlet cock on the receiver.

2. Compressor does not start, makes a buzzing sound followed by a clicking noise (cannot start against high pressure):

- a) Leaky non-return valve. Remove the flexible pressure pipe and clean to find out whether air leaks from the valve. If so, clean or replace.

3. Vacuum pump operates, but does not evacuate sufficiently:

- a) Check the piston gaskets. Replace, if necessary.
- b) Defective valve plate. Contact your JUN-AIR distributor.
- c) Inlet cock defective. Replace.

If the vacuum pump is fitted with receiver, vacuum filter or vacuum gauge:

- d) Vacuum filter is clogged. Clean or replace filter element.
- e) Leaks in fittings, hoses or equipment. May be located by closing inlet cock. Dead-end vacuum is to correspond to the charts on page 41.
- f) Vacuum gauge defective. Replace.

4. Compressor works, but pressure does not increase:

- ⊖ a) Intake filter clogged. Replace.
- b) Leaks in fittings, tubes or pneumatic equipment. Check with soapy water or by letting unit stay overnight with disconnected mains. Pressure drop is not to exceed 1 bar.
- c) Check the piston gaskets. Replace, if necessary.
- d) Defective valve plate. Contact your JUN-AIR distributor.
- e) Failure in non-return valve which is creating a flow restriction.
- ⊖ f) Leaks in fittings, tubes or pneumatic equipment. Check unit by letting unit stay overnight with disconnected mains. Vacuum switch is not to drop to 0 bar.

5. Loud noise from compressor:

- ⊖ a) Dirt or failure in non-return valve. Clean or replace.

6. Motor gets very hot:

- a) Leaks. See point 4b.
- b) Clogged vacuum filter. See point 4a.
- c) Too high ambient temperature. Ensure adequate ventilation if the vacuum pump is installed in a cabinet.
- d) Overloaded. Ensure compressor is correct model for work load.

7. Compressor starts when no air is being tapped:

- ⊖ a) Leaks. See point 4b.

8. Vacuum pump starts when no vacuum is being used:

- ⊖ a) Leaks. See point 4f.

9. Vacuum pump starts and stops more frequently than usual:

- ⊖ a) Leaks. See point 4f.

10. Vacuum pump does not start at minimum vacuum or does not stop at stop vacuum:

- ⊖ The following points only apply to units with vacuum gauge:
- a) Vacuum switch defective. Replace.
- b) Vacuum switch is not adjusted correctly. Adjust correctly.

11. Compressor does not switch on against pressure or does not switch off at max. pressure:

- ⊖ a) Defective pressure switch. Replace.

Test the non-return valve

- Once a year.
- Switch off compressor on the main switch and pull out the plug.
- Pull the ring at the end of the safety valve (fig. 17).

Warning:

Loud noise!

- Let pressure fall to 0 bar. The pressure is to be read on the pressure gauge for receiver pressure (fig. 12).
- If the compressor is mounted with a TÜV-approved safety valve, the receiver is emptied by loosening the screw at the end of the safety valve (fig. 19).
- Dismount non-return valve from receiver (fig. 25).
- Disassemble non-return valve and remove O-ring from piston (fig. 26 and 27).
- Clean non-return valve.
- Mount new O-rings and re-assemble non-return valve (fig. 28). (Two O-rings in non-return valve)

Preventive maintenance (all)

	Weekly	Monthly	Annually
Drain condensate from air receiver. If equipped with autodrain, this will take place automatically, however, drain bottle has to be emptied.	•		
If compressor is fitted with outlet filter, check and empty for water by pressing the black button at the bottom. If fitted with autodrain, this will take place automatically.	•		
Check motor, air tubes, hoses and equipment for leaks, and check the pumping time.		•	
Clean unit or wipe with a soft, damp cloth. If necessary, use paraffin on rag to remove sticky adhesions. Dust and dirt prevent cooling.		•	
Check intake filter. Clean or replace if necessary. Min. annually or every 2000 hours whichever comes first.			•
Replace silencer (intake filter), if necessary (fig 13 - Loosen air filter using e.g. a large coin (1)) + (14 - Replace and tighten new air filter using a coin (2)).			•
Check the O-ring in the non-return valve and replace, if necessary. Note! Empty receiver of air before dismounting.			•
Check filter and filter elements for optimum efficiency.			•
Test the safety valve by gently pulling the ring with pressure in the receiver.			•

Replace coal brushes on DC motors 2000 hours

Please note: Vacuum pumps without air receiver are not supplied with vacuum filter.

Check the pumping time (Comp.)

The pumping time indicates the condition of the compressor provided that there are no leaks in the system where the compressed air may leak. Test the compressor as follows:

1. Empty the air receiver of compressed air (the pressure gauge shows 0 bar).
2. Close the outlet on the air receiver and check that the drain cock is closed.
3. Start the compressor and note how long it takes until it switches off.

Ensure that the pressure in the air receiver is 8 bar/120 psi as deviations may indicate the wrong results (see Technical Specifications for pumping time).

Important!

Always test the compressor when cold as the time indicated refers to the pumping time of a cold compressor. The pumping time of a warm compressor is much longer and consequently, the result would be misleading.

Pressure vessel

Pressure tested at:	4-25 litre:	24 bar
	40-50 litre:	18.3 bar

Directions for use

Application	Receiver for compressed air.
Receiver specifications	See name plate.
Installation	Tubes, etc. must be installed with suitable materials.
Placement	Observe the working temperature of the receiver. Ensure sufficient room for inspection and maintenance. The receiver must be kept in a horizontal position.
Corrosion protection	The surface treatment must be maintained as required. Internal inspection at least every 5 years. Drain condensate at least once a week.
Alternation/repair	No welding must be made on pressurised parts.
Safety valve	Ensures that PS will not be exceeded. Never adjust to a higher pressure than PS. The capacity of the valve must be calculated in accordance with the volume of air supplied by the compressor. (PS = Maximum working pressure of the receiver)

Declaration of Conformity

NOTE: The declaration of conformity is only valid for units operating at 230 V/50 Hz, 3x400 V/50 Hz, 12 V DC or 24 V DC.

The manufacturer, JUN-AIR International A/S, declares that the products mentioned in this manual are in conformity with:

- 87/404/EEC - 90/488/EEC - 93/68/EEC Council Directive relating to Simple Pressure Vessels
- 98/37/EC Machinery Directive
- 89/336/EEC EMC Directive
- 73/23/EEC Low Voltage Directive

Flemming Petersen
Quality Manager

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* Neutral is required

Technical modifications reserved

Model		Motor OF301 vacuum										
Voltage	Volt	12	24	120	230	230	3x400*	3x400*				
Frequency	Hz	DC	DC	60	50	60	50	60				
Power	HP	0,33	0,33	0,33	0,33	0,33	0,33	0,33				
	kW	0,24	0,24	0,24	0,24	0,24	0,24	0,24				
Displacement	l/min	69	69	68	54	68	54	68				
	CFM	2,44	2,44	2,40	1,91	2,40	1,91	2,40				
Dead end	l/min	100	100	100	100	100	100	100				
	CFM	3,53	3,53	3,53	3,53	3,53	3,53	3,53				
Max. current	Amps	27	14	4,2	2,1	2,2	0,8	0,8				
Weight	Kg	10	10	10	10	10	10	10				
	Lbs	22	22	22	22	22	22	22				
Dimensions (l x w x h)	mm	281 x 139 x 239										
	Inch	11,1 x 5,5 x 9,4										
Noise level @ 1 m	dB(a)	63	63	61	60	61	60	61				

Model		Motor OF302 vacuum										
Voltage	Volt	120	230	230	3x400*	3x400*						
Frequency	Hz	60	50	60	50	60						
Power	HP	0,60	0,60	0,60	0,60	0,60						
	kW	0,44	0,44	0,44	0,44	0,44						
Displacement	l/min	138	108	138	108	138						
	CFM	4,87	3,81	4,87	3,81	4,87						
Dead end	l/min	100	100	100	100	100						
	CFM	3,53	3,53	3,53	3,53	3,53						
Max. current	Amps	6,6	3,4	3,8	1,4	1,4						
Weight	Kg	13	13	13	13	13						
	Lbs	29	29	29	29	29						
Dimensions (l x w x h)	mm	312 x 139 x 238										
	Inch	12,3 x 5,5 x 9,4										
Noise level @ 1 m	dB(a)	61	60	61	60	61						

Motor size		Motor OF322										
Voltage	Volt	120	230	230	3x230*	3x230*	3x400*	3x400*				
Frequency	Hz	60	50	60	50	60	50	60				
Power	HP	0,60	0,60	0,60	0,60	0,60	0,60	0,60				
	kW	0,44	0,44	0,44	0,44	0,44	0,44	0,44				
Displacement	l/min	84	67	84	67	84	67	84				
	CFM	2,97	2,37	2,97	2,37	2,97	2,37	2,97				
FAD @ 8 Bar ***	l/min	44	40	44	40	44	40	44				
	CFM	1,55	1,41	1,55	1,41	1,55	1,41	1,55				
Max. pressure**	Bar	12	12	12	12	12	12	12				
	PSI	175	175	175	175	175	175	175				
Max. current	Amps	8,4	3,8	4,4	3,0	2,5	1,6	1,4				
Weight	Kg	13	13	13	13	13	13	13				
	Lbs	29	29	29	29	29	29	29				
Dimensions (l x w x h)	mm	320 x 140 x 240										
	Inch	12.5 x 5.5 x 9.5										
Noise level @ 1 m	dB(a)	72	72	72	72	72	72	72				

* Neutral is required

Technical modifications reserved

Motor size		Motor OF311			Motor OF312							
Voltage	Volt	120	230	230	120	230	230					
Frequency	Hz	60	50	60	60	50	60					
Power	HP	0,33	0,33	0,33	0,60	0,60	0,60					
	kW	0,24	0,24	0,24	0,44	0,44	0,44					
Displacement	l/min	84	67	84	168	134	168					
	CFM	2,97	2,37	2,97	5,93	4,73	5,93					
FAD @ 8 Bar ***	l/min	34	26	34	64	58	64					
	CFM	1,20	0,92	1,20	2,26	2,05	2,26					
Max. pressure**	Bar	6	6	6	6	6	6					
	PSI	90	90	90	90	90	90					
Max. current	Amps	4,7	2,6	2,7	8,0	3,4	3,8					
Weight	Kg	10	10	10	13	13	13					
	Lbs	22	22	22	29	29	29					
Dimensions (l x w x h)	mm	230 x 140 x 240			320 x 140 x 240							
	Inch	9.0 x 5.5 x 9.5			12.5 x 5.5 x 9.5							







Model			OF301 - 4B	OF302 - 4B	OF302 - 15B	
Tank size	liter		4	4	15	
	US gallon		1,1	1,1	4,0	
Weight	kg		19	22	25	
	lbs		42	49	55	
Dimensions (l x w x h)	mm		390 x 310 x 350	390 x 320 x 350	380 x 380 x 530	
	Inch		15,5 x 12,5 x 13,5	15,5 x 12,5 x 13,5	15,0 x 15,0 x 21,0	
Pumping time (0-8 bar/0-120 psi)	@ 50 Hz	sec.	90	45	165	
	@ 60 Hz	sec.	75	40	135	
Noise level @ 1 m	@ 50 Hz	dB(a)	65	65	65	
	@ 60 Hz	dB(a)	66	66	66	




Model			OF302 - 25B	OF302 - 25BD2**	2*OF302 - 40B	2*OF302-40BD2**
Tank size	liter		25	25	40	40
	US gallon		6,6	6,6	10,6	10,6
Weight	kg		28	33	46	59
	lbs		62	73	101	130
Dimensions (l x w x h)	mm		380 x 380 x 610	510 x 460 x 610	560 x 450 x 610	560 x 540 x 630
	Inch		15,0 x 15,0 x 24,0	20,0 x 18,0 x 24,0	22,0 x 18,0 x 24,0	22,0 x 21,0 x 25,0
Pumping time (0-8 bar/0-120 psi)	@ 50 Hz	sec.	270	325	220	260
	@ 60 Hz	sec.	225	250	180	220
Noise level @ 1 m	@ 50 Hz	dB(a)	65	65	68	68
	@ 60 Hz	dB(a)	66	66	69	69

** Displacement is reduced by approx. 18-20% on units with dryer (D). Min. pressure required to operate dryer: 6 bar
 Technical modifications reserved.


The figures in the tables are based on the unit working in fairly clean air, at ambient temperatures around 20-40°C, relative humidity near 50% and operating at sea level. The performance of the product will be adversely effected at high altitudes. For further information on the wirings isolations class: See data table on the motor.

Translations					
English	German	French	Spanish	Dutch	Danish
Voltage	Spannung	Voltage	Voltaje	Voltage	Spænding
Frequency	Frequenz	Fréquence	Frecuencia	Frequentie	Frekvens
Motor HP	Motor HP	Moteur CV	Motor CV	Motor HP	Motor HK
Displacement	Ansaugleistung	Débit	Aire aspirado	Capaciteit	Ydelse
Max. pressure	Max. Druck	Pression de service max.	Presión de régimen máx.	Max. druk	Max. arbejdstryk
Max. current	Stromverbrauch	Consommation	Corriente máxima	Max. stroom	Strømforbrug
Tank size	Behältervolumen	Volume réservoir	Volumen de tanque	Tankvolume	Beholderstørrelse
Weight	Gewicht	Poids	Peso	Gewicht	Vægt
Dimensions (l x w x h)	Abmessungen (l x b x h)	Dimensions (l x p x h)	Dimensiones (l x a x h)	Afmetingen (l x w x h)	Dimensioner (l x b x h)
Noise level	Schallemissionen	Niveau sonore	Nivel de ruido	Geluidsniveau	Lydniveau
Pumping time	Pumpzeit	Temps de refoulement	Tiempo de bombeo	Pomptijd	Oppumpningstid
Neutral is required	Null-leiter ist erforderlich	Neutre nécessaire	Neutro necesario	Neutraal noodzakelijk	N-leder kræves
Displacement is reduced by approx. 18-20% on units with dryer (D). Min. pressure required to operate dryer: 6 bar	Bei Kompressoren mit Adsorptionstrockner reduziert sich die effektive Luftliefermenge um 18-20% (D). Mindestarbeitsdruck für den Luft-trockner beträgt 6 bar	Le débit est réduit de 18-20% pour les unités avec sécheur d'air (D) Pression min. 6 bar	Le capacidad se reduce con 18-20% para las unidades con secador de aire (D)	Bij systemen met droger is de capaciteit ca. 18-20% lager (D). Min. benodigde druk voor de droger: 6 bar	Ydelsen reduceres med ca. 18-20% på kompressorer med tørrer (D). Min. tryk til drift af tørrer: 6 bar
Technical modifications reserved	Technische Änderungen vorbehalten	Droits réservés pour modifications techniques	Reservamos el derecho a cambiar estas especificaciones técnicas sin previo aviso	Technische wijzigingen voorbehouden	Ret til ændringer forbeholdes


Dimensions - OF301 12/24V DC


0070011

Dimensions - OF301 12/24V DC, horizontally mounted


0070015

Dimensions - OF301 50Hz



0070013

Dimensions - OF301 60Hz



0070010

Dimensions - OF301 50Hz, horizontally mounted


0070012

Dimensions - OF301 60Hz, horizontally mounted


0070014

Dimensions - OF302


0080010

Dimensions - OF302, horizontally mounted


0080011

Electrical diagrams - OF301/302 - 3x400V / 50/60Hz


6072864

Electrical diagrams - OF301 - 12/24V DC with motor protection


6072866

Electrical diagrams - OF301 - 12/24V DC without motor protection

WARNING: NO MOTOR PROTECTION

6072867


For internal motor wiring, please refer to operating manual

SPECIFICATION


Motorspec.: OF301V, 100V~50Hz - 3.8A
OF301V, 100V~60Hz - 3.8A
OF301V, 120V~60Hz - 3.1A
OF301V, 200V~50Hz - 1.9A
OF301V, 200V~60Hz - 1.8A
OF301V, 230V~50Hz - 1.5A
OF301V, 240V~60Hz - 1.8A

Motorspec.: OF302V, 100V~50Hz - 5.5A
OF302V, 100V~60Hz - 6.4A
OF302V, 120V~60Hz - 6.6A
OF302V, 200V~50Hz - 2.9A
OF302V, 200V~60Hz - 3.7A
OF302V, 230V~50Hz - 2.9A
OF302V, 240V~60Hz - 3.2A


-  :Earth
-  : Terminal clip
-  :Capacitor
-  :Motor protection
-  :Box

For internal motor wiring, please refer to operating manual

For internal motor wiring, please refer to operating manual

For internal motor wiring, please refer to operating manual.


Spare parts

2881100	Motor compressor OF301 230V/50Hz
2881200	Motor OF301 120V/60Hz
2881206	Motor OF301 230V/60Hz
2881302	Motor OF301 3x400V/50Hz without motor protection
2882100	Motor OF302 230V/50-60Hz
2882200	Motor OF302 120V/60Hz
2882400	Motor compressor OF302 3x400V/50-60Hz
3050099	Receiver 4l with inspection plug CE, silver, internally coated
3210093	Receiver 15l w/inspection plug and internally coated
3210097	Receiver 15l with inspection plug and internally coated, RAL9003, ASME
3410095	Receiver 25l with inspection plug, internally coated, ASME, silver
3410098	Receiver 25l CE silver with inspection plug and internally coated
3512098	Receiver 40l silver w/inspection plug internally coated CE/ASME
4071020	Regulator R07-280-RNMG, 10 bar
4071030	Filter regulator 5um w/manual drain
4071125	Regulator 16 bar w/gauge
4082524	Dryer f/Basic OF302-25BD2 - 8 bar - type OX (size 2)
4082928	Dryer f/2xOF302-40BD2/OF1201-40BD2 - OX - 8 bar
4146458	Rapid fitting Ø8mmxØ8mm elbow
5030000	Pressure switch MDR 2/11 compl.
5075010	Pressure switch MDR 21/11 UL compl.
5130000	Gauge Ø40 0-16 bar 1/8" down
5235805	Capacitor, run complete 10uF/230V AC f/OF300 (30 cm)
5237505	Capacitor, run complete 25uF/230V AC f/OF300 (35 cm)
5237506	Capacitor, run complete 25uF/230V AC f/OF300 (60 cm)
5240605	Capacitor, run complete 60uF/120V AC f/OF300 (35 cm)
5320226	Handle, multi complete
5340001	Handle f/6-15/6-25 complete
5412400	Intake filter f/OF301 and OF302 motors
5412401	Intake filter f/OF300 motors - short
5414500	Non return valve
5416100	Safety valve 12 bar / 177 psi
5417000	Safety valve 9 bar / 135 psi
5417700	Drain cock 1/4"
5418000	Drain cock 1/4" 4l
5419500	Drain cock 1/4" 15l
5420000	Drain cock 1/4" 25l
5421001	Drain cock 1/4" 40l (new)
5425700	Safety valve TÜV 12.2 bar
5427000	Connecting piece 1/4"

5445200	Radiator & fan M2000 with cam timer 230V CE/UL
5445300	Radiator & fan w/camtimer M2000 120V CE/UL
5459500	End cover f/OF301 and OF302 motors
5470500	Spare parts kit f/replacement of valve plate f/OF301 and OF302 motors
5470501	Spare parts kit f/replacement of valve plate f/OF301V and OF302V
5470510	Spare parts kit f/replacement of valve plate f/OF311 and OF312 motors Ø70
5470520	Spare parts kit f/replacement of valve plate f/OF301/OF302 (ø63)
5470521	Spare parts kit f/replacement of valve plate f/OF301V/OF302V
5470522	Spare parts kit f/replacement of valve plate f/OF312/OF322 low pressure (ø70)
5470523	Spare parts kit f/replacement of valve plate f/OF322 motors (high pressure) Ø63
5499800	Connecting rod/cylinder complet f/OF301/OF302
5499801	Connecting rod/cylinder complet f/OF301V/OF302V
5499810	Connecting rod/cylinder complet f/OF311/OF312/OF322 (ø70)
5499811	Connecting rod/cylinder complet f/OF322/332 (Ø63)
5612300	Drain bottle 1 l multi supplied separately incl. bracket for pipe mounting and magnet
5960000	Cable relief
5986000	Fan blade f/OF301/OF302 black CW
5986010	Fan blade f/OF301/OF302 grey CCW
6236952	Gasket, valve plate/cylinder head OF300 (silikone)
6241800	O-ring f/1" plug - 32x5 mm
6243710	O-ring f/plug Ø11.2x2.4 f/OF300
6245800	O-ring Ø47x5 f/1 1/2" plug
6252010	Socket f/handle 3-1.5/6-4 Ø22
6253000	Rubber base f/handle Ø25x30 M8x8
6253200	Rubber base Ø50
6253300	Rubber base complete f/3-4000 models
6253701	Vibration damper f/OF300 / 2-600 Ø20x30
6254300	Vibration damper f/OF301 230V
6290200	Plastic handle black Ø22 mm
6317000	Nut M12 DIN934-8
6340500	Unbraco plug 1/8"
6357800	Washer Ø28.5xØ13x1.5 FZB
6392155	Screw torx M4x20 PH
6420000	Flex hose 1/8" 20 (22) cm
6433000	Flex hose 1/8" 33 (35) cm
6451016	Teflon hose 1/4" 16 cm
6451028	Teflon hose 1/4" 28 cm
6452024	Teflon hose 1/4" 24 cm with elbow
6452028	Teflon hose 1/4" 28 cm with elbow
6452042	Teflon hose 1/4" 42 cm with elbow
6452065	Teflon hose 1/4" 65 cm with elbow


6470320	Unloader complete 20 cm with noise reduction
6470326	Unloader complete 26 cm with noise reduction
6470333	Unloader complete 33 cm with noise reduction
6482000	Flex hose 1/4" 23 (25) cm
7020000	Double nipple 1/8"
7024000	Double nipple 1/4" L= 26 mm
7024002	Double nipple 1/4" w/o-ring track
7025000	Double nipple 1/4"x1/8"
7070000	Connection piece f/non-return valve
7070600	Connection piece f/non-return valve model 12-25
7071000	T-piece 1/8"
7071400	T-piece Condor
7071500	T-piece 1/4" (3 x internal)
7156000	Extension piece f/cock short
7156100	Extension piece f/cock long
7157001	Extension piece w/through 1/4"
7164200	Cross connector
7164900	Plug 1" f/inspection
7166033	Plug 1 1/2" f/inspection ASME
7172800	Spacer f/600-40/50B
7180000	Elbow 1/4" int./ext. KRG
7181002	Elbow 1/4" x 1/4" f/outlet
7190000	Elbow 1/8"
7565000	Bushing 1/4" ext. x 1/8" int.
7566000	Bushing 1/4" x 1/4"
7566200	Bushing 1/4" x 1/4" L = 35 mm
8111511	Motor foundation f/12-40, painted

Spare parts OF301 12/24V DC


0070002

Spare parts OF301


0070000

Spare parts OF311


0000016

Spare parts OF302


0080000

Spare parts OF312


0080014

Spare parts OF322


0080015

Sparte parts model OF301-4B


0025301

Sparte parts model OF302-4B


0025302

Spare parts model OF302-15B


Spare parts model OF302-25B


0025402

Spare parts model OF302-25BD2


0025403

Spare parts model 2xOF302-40B

Spare parts model 2xOF302-40BD2

0025405

Fig. 1a


Fig. 1b

Fig. 1c

Fig. 2


Fig. 3


Fig. 4a

Fig. 4b

Fig. 5

Fig. 6

Fig. 7

Fig. 8


Fig. 9


Fig. 10


Fig. 11


Fig. 12


Fig. 13

Fig. 14

Fig. 17


Fig. 25


Fig. 26

Fig. 27

Fig. 28

JUN-AIR®

P.O. Box 97
Benton Harbor,
Michigan 49023-0097
USA

www.jun-air.com

Phone: 269-934-1216
Fax: 269-927-5725
E-mail: jun-air@idexcorp.com

GAST

A UNIT OF IDEX CORPORATION

IDEX
IDEX CORPORATION